

Characterizing immunotherapy-induced lymphocyte infiltration at the single patient level using CANscript™, an ex-vivo human tumor model

Munisha Smalley¹, Basavaraja U Shanthappa¹, Hans Gertje¹, Mark Lawson¹, Baraneedharan Ulaganathan¹, Allen Thayakumar¹, Laura Maciejko¹, D.C. Doval², Anurag Mehta², Padhma Radhakrishnan¹, Pradip Majumder¹, Aaron Goldman¹






1. Integrative Immuno-Oncology Center, Mitra Biotech. Woburn, MA 01801, USA 2. Rajiv Gandhi Cancer Institute, New Delhi, India.

Background: The presence and activity of lymphocytes within the tumor is critical for clinical response to cancer immunotherapy, such as immune checkpoint blockade. Tumors with poor T-cell inflamed phenotypes, often referred to as a 'cold' tumor, is associated with modest clinical response. High baseline infiltration of effector T-cell lymphocytes is considered 'hot', and patients are predicted to respond more favorably to treatment. However, patient-to-patient response and durability remains highly variable. There is an urgent gap in available methods to study lymphocyte infiltration, trafficking and spatial heterogeneity induced by different cancer immunotherapies in individual patients. Moreover, there is a poor correlation between therapy-induced lymphocyte infiltration with clinical response, which could be shaped using personalized approaches to therapy.

Methods: Here, we used CANscript™, an ex-vivo human tumor model that recapitulates and preserves the native, patient-autologous tumor microenvironment, including autologous patient-derived peripheral blood mononucleated cells (PBMC). Utilizing tissue from breast cancer patients classified as either 'cold' (N=5) or 'hot' (N=5), we studied lymphocyte infiltration under pressure of α -PD-1 immune checkpoint blockade (pembrolizumab) over a 72h time course compared to standard of care (SOC) drugs such as docetaxel. Using fluorescent labelling and flow cytometric analysis we characterized infiltrating lymphocytes, studying the role of T-cell repertoires under different environmental and immunotherapy pressures. We coupled these analyses with multiplex immunohistochemistry (CD3⁺, CD4⁺, CD8⁺) to map spatial heterogeneity of tumor cells and lymphocytes before and after treatment, ex-vivo.

Results: We determined that immune checkpoint blockade induced unique patterns of migration and infiltration of effector T-cells (T_{eff}) and T-regulatory (T_{reg}) cells in 'hot' vs 'cold' tumors. Furthermore, we determined that, in some instances, 'cold' tumors can be driven towards a 'hot' phenotype characterized by trafficking of active immune lymphocytes following treatment, which corresponded to differential ratio of T_{eff} to T_{reg} compared to baseline.

Concluding remarks: Taken together, these data demonstrate the utility of CANscript™ as a platform to characterize response to immunotherapy in a spatial context, providing insight into the migratory patterns of immune cell subsets at the individual patient level. Such an advance in our preclinical methods to study immuno-modulators may help guide treatment decisions for clinicians while simultaneously functioning as a platform to study and discover mechanisms of clinical efficacy for emerging drug combinations.



Conclusions: Dynamic, spontaneous infiltration of lymphocytes into the tumor (aka TILs) are a key factor for immune-based tumor rejection. However, little is known about the role of immunotherapy or other anticancer agents as they modulate TILs. CANscript™ enables a platform to study the antitumor effect of different therapies along with stochastic lymphocyte infiltration at the individual patient-level. By integrating our algorithm-driven predictive clinical assessment (M-Score) we are able to provide unique understanding towards the role of dynamic TIL with clinical response to therapy. Our data demonstrate that high inter-patient variability in immune modulation, and response to immuno-modulators, suggests high importance of integrating a multi-dimensional approach to personalized cancer medicine.

