


Introduction

Evaluating genotoxicity is a crucial part of chemical risk assessment and includes tests for the induction of gene mutations, chromosomal aberrations, and numerical chromosomal changes. Traditionally, mutagenicity is measured as a result of phenotypic changes but this underestimates the true mutation frequency (MF). A sequencing approach facilitates the detection of all mutations but historically suffered from poor precision due to inherent amplification errors when attempting to quantify $< 1 \times 10^{-7}$ mutations/kb.


However, with the introduction of error-corrected sequencing (ecSeq) methods, it is now possible to resolve ultra-low MFs. Single-molecule mutation sequencing (SMM-seq) is a highly sensitive technique for detecting single nucleotide variants (SNVs). SMM-seq utilizes Rolling Circle Amplification (RCA), which amplifies linked strands of each DNA fragment into a concatenated single-stranded DNA product. These ssDNA contigs act as a template for the sequencing library. This allows for the direct comparison of multiple copies of replicated strands, reliably distinguishing genuine mutations from sequencing artifacts.

ToxTracker is a new approach methodology (NAM) that evaluates the induction of pathway-specific GFP-reporter genes to discriminate between primary and secondary genotoxicity, thereby distinguishing the likelihood of DNA reactivity. Specifically, the assay identifies secondary genotoxicity caused by oxidative stress or protein damage and can distinguish aneugens from clastogens. Here, we combined ToxTracker with SMM-seq to evaluate the genotoxic and mutagenic mode of action of nine substances.


The ToxTracker® assay



Single-molecule mutation Sequencing



ToxTracker predicts the genotoxic mode of action



Mutational signatures reveal the mode of action of genotoxic compounds

