



Optimisation of *in vitro* metabolism using S9 liver extract in ToxTracker®

R Derr¹, N Moelijker¹, Lorrie Boisvert², I Brandsma¹, Paul White², G Hendriks¹

1. Toxys B.V., Leiden, the Netherlands
2. Health Canada, Ottawa, Canada


Introduction

Progenotoxic agents, such as Aflatoxin B1 and Benzo[a]pyrene, only become DNA reactive after metabolism by detoxification enzymes. *In vivo* metabolism takes place in for example the liver, bone marrow or the lungs. In *in vitro* genotoxicity assays such as ToxTracker, metabolism can be included by using S9 liver extract from rat or hamster. To induce the expression of detoxification enzymes, animals were treated with Aroclor-1254 or a mixture of phenobarbital and 5,6-benzoflavone.


In ToxTracker, we tested two different rat S9 liver extracts and compared their capacity to metabolise 20 progenotoxic compounds. Furthermore, we optimised the S9 concentration that was used in the assay to allow longer exposures in the presence of S9 without inducing cytotoxicity. The potency of the different S9 extracts and impact of the modifications of S9 concentration and exposure times were analysed by comparing the LOAEL.

The ToxTracker reporter assay


- ToxTracker is a stem cell-based reporter assay consisting of 6 GFP reporter cell lines.
- GFP-tagged biomarkers are activated upon specific cellular responses to DNA damage or other stress.
- Activation of the biomarkers is analysed using flow cytometry.



S9 toxicity



Optimisation of S9-mediated drug metabolism



Summary S9 metabolism efficiency

Table 1: Lowest observed adverse effect concentration to induce genotoxicity or oxidative stress in the absence and presence of S9.

| Compound | -S9 | | 3h. 1% S9 | | 24h. 0.25% S9 | |
|---------------------------------------|-------|-------|-----------|-------|---------------|-------|
| | GTX | OX | GTX | OX | GTX | OX |
| Aflatoxin B1 (AFB1) | - | - | 0.63 | 2.5 | 0.63 | 1.25 |
| Benzo[a]pyrene (B[a]P) | - | - | 12.5 | 12.5 | 6.25 | 6.25 |
| Cyclophosphamide (CPA) | - | - | 6.25 | 12.5 | 12.5 | 12.5 |
| PhiP | - | - | 3.91 | 7.81 | 0.98 | 1.95 |
| Tryptophan-P-2 | - | - | 2.44 | 4.88 | 1.22 | 1.22 |
| 7,12-Dimethylbenz(a)anthracene (DMBA) | - | 15.63 | 1.95 | 3.91 | 3.91 | 3.91 |
| 2-aminoanthracene (2-AA) | - | 9.77 | 19.53 | 19.53 | 9.77 | 9.77 |
| IQ | 1000 | - | 62.5 | 250 | 62.5 | 62.5 |
| 3-methylcholanthrene (3-MC) | 625 | 250 | 15.62 | 31.25 | 31.25 | 62.5 |
| MeIQ | 500 | - | 62.5 | 62.5 | 62.5 | 500 |
| 2-Acetylaminofluorene (2-AAF) | 100 | 100 | 125 | - | 100 | 50 |
| 1,2-Diphenylhydrazine (1,2-DPH) | 312.5 | 312.5 | - | 39.06 | 312.5 | 39.06 |
| Disperse Orange | 31.25 | - | 62.5 | - | - | - |
| 1,3-diphenyltriazine (1,3-DPT) | 39.06 | 39.06 | 39.06 | 78.13 | 39.06 | 39.06 |
| Phenolphthalein | 39.06 | 39.06 | 312.5 | 312.5 | 39.06 | 19.53 |
| Sodium arsenite | - | 125 | - | - | - | 1.25 |
| ortho-aminoazotoluene (OAT) | - | 78.12 | - | 312.5 | - | 78.13 |
| Dimethyl yellow | - | 31.25 | - | - | - | 15.63 |
| Isoprene | - | - | - | - | - | - |
| Hexamethylphosphoramide (HMPA) | - | - | - | - | - | - |

All compounds were tested in the ToxTracker assay in the absence of S9, in the presence of 1% S9 for 3h, or 0.25% S9 for 24h. The LOAEL was determined based on ToxTracker data to compare the metabolic activation of the substances. For genotoxicity assessment, the LOAEL of Rtkn-GFP and Bcl2-GFP were determined. For the oxidative stress LOAEL, activation of the Srxn1-GFP and Blvr-GFP reporters were assessed. The LOAEL was defined as the concentration at which a 2-fold induction of either reporter was observed.

Aroclor-1254 vs. Phenobarbital induced S9




Figure 3: Differential activation of the ToxTracker reporters in presence of 0.25% aroclor-1254 of Phenobarbital induced S9 metabolising system.

| Compound | Aroclor-1254 S9 | Phenobarbital S9 |
|---------------------------------------|-----------------|------------------|
| Aflatoxin B1 (AFB1) | 0.625 | 1.25 |
| Benzo[a]pyrene (B[a]P) | 6.25 | 6.25 |
| Cyclophosphamide (CPA) | 12.5 | 25 |
| PhiP | 0.98 | 1.95 |
| Tryptophan-P-2 | 1.22 | 1.22 |
| 7,12-Dimethylbenz(a)anthracene (DMBA) | 3.90 | 15.62 |
| IQ | 62.5 | 500 |
| 3-methylcholanthrene (3-MC) | 31.25 | 62.5 |
| MeIQ | 62.5 | 500 |
| 2-Acetylaminofluorene (2-AAF) | 100 | 100 |
| 1,2-Diphenylhydrazine (1,2-DPH) | 312.5 | 312.5 |
| 1,3-diphenyltriazine (1,3-DPT) | 39.06 | 39.06 |
| Phenolphthalein | 39.06 | 78.13 |
| 2-aminoanthracene (2-AA) | 9.77 | - |
| Benz[a]anthracene | 62.5 | - |