


TubulinTracker, a Novel In Vitro Reporter Assay to Study Intracellular **Microtubule Dynamics and Cell Cycle Progression**

Introduction


Aneuploidy, the presence of an abnormal number of chromosomes, can be caused by any process that interferes with chromosome segregation during mitosis, including microtubule disruption or inhibition of the Aurora A/B/C cell cycle kinases. To detect aneugenicity as a result of chemical exposure, typically the micronucleus assay is used. However, both broken DNA fragments caused by clastogenic agents and mis-segregated chromosomes caused by aneugenic agents can lead to the formation of micronuclei.


Here, we established a reporter assay to study the effect of substances on tubulin stability, to provide more insight into the cause of aneuploidy. The use of a GFP-tagged tubulin reporter allows the direct visualisation of microtubuli, using either live cell microscopy to follow the dynamics during the cell cycle or flow cytometry.

To measure the effect of substances on microtubule stability, cells were treated, fixed and the GFP signal was quantified using flow cytometry. In order to assess the effect of the compounds on the cell cycle distribution, a DNA staining (Hoechst) was included. The mechanistic insight into the MOA provided by this assay is important for hazard identification and part of weight of evidence approaches.





I. Brandsma, N. Moelijker, M. Hoogenboom-Geijer, R.Derr and G. Hendriks






decreased. After both treatments, the number of mitotic cells stronly increased. White arrows indicate mitotic cells.





Toxys B.V., Leiden, The Netherlands

percentage of cells at each cell cycle phase, after 4h of treatment with the indicated compounds. <u>Right</u> panels, the analysis of GFP-tubulin signal in G2/M cells only. Exposure to all substances lead to an increase in G2/M phase cells. Tubulin stabilisers significantly increased, while tubulin destabilisers significantly decreased the GFP-Tubulin signal. Bars represent the mean for 3 biological replicates, error bars show SEM. *p<0.05 **p<0.005 treated vs DMSO control (t-test).



Overview of Tubulin stabilising and destabilising compounds. The mES GFP-Tubulin assay was validated using 23 compounds: 3 microtubule stabilisers, 8 destabilisers and 12 agents not known to affect microtubuli. GFP-Tubulin cells were treated for 4h, using 10 concentrations of the test substances. To classify agents, the area under the curve (AUC) was calculated for the exposures, normalising both the exposure concentration as well as the fold-change in GFP-signal. Compounds stabilising microtubuli showed a large increase in AUC, while tubulin destabilising compounds showed a decrease in AUC. For compounds not affecting tubulin stability, the AUC was around 1 (dotted line).

Conclusions

• Microtubule dynamics were studied in the context of a living cell, using physiological expression levels of the GFP-Tubulin.

Using a stable GFP-Tubulin reporter cell line microtubule poisons can be reliably identified.

- Exposure to other agents affecting mitotic progression or causing other types of damage does not affect the GFP-Tubulin signal.
- The assay provides insight into the MOA of aneugenic agents